



Monitoring your Blood Glucose levels with your GM3+ Glucometer (*Blood Glucose Meter*) can help you take better care of your diabetes. It will help you learn how food, activity levels, stress, medicine and insulin change your blood sugar level. Just a tiny droplet of Blood and within 20 seconds you will have your Blood Glucose values.

**Make your Smart Choice Today !**



# easier than ever before

# 1



Insert the strip into the GM3+ monitor and switch ON the monitor.

# 2



Apply a tiny sample of blood at the test area of the test strip.

# 3



In 20 seconds, your GM3+ Monitor displays the Blood Glucose values in mg/dl.



**BIG** FEATURES  
SMALL PRICE

Designed  
& Made  
in India

Small  
&  
Portable

Large  
Display  
Characters

Test Results  
in  
20 seconds

Wide  
Test Range  
20 mg/dl  
to 600 mg/dl

Large 75  
Results  
Memory

Extended  
Battery Life

User  
Replaceable  
AAA Battery

Strips  
Specially  
Packed for  
extended life

Memory  
Smart Tags  
for Before  
& After Meal

Displays  
Average  
of ALL  
Readings

Displays  
Before  
& After Meal  
Averages

Support by  
Email, Fax  
Telephone  
or Online

Buy Online  
or  
Order by  
Telephone

Buy from  
the extensive  
Dealers  
Network

## SPECIFICATIONS:

Product:	GM3+ Blood Glucose Monitor
Application:	To determine blood glucose in whole blood by amperometric method
Size:	120mm x 70mm x 24mm
Weight:	110gms (including batteries)
Operating Conditions:	Temperature: 5°C to +45°C Humidity: Less than 98%
Power Supply:	User Replaceable Two 1.5-Volt Alkaline AAA batteries(Type LR 03)
Battery Life:	Around 10000 tests (2 minutes per test).
Display:	Special Large Size Display with Blood Glucose display characters 12 mm in height.
Measuring range:	20 mg/dl to 600 mg/dl
Measuring time:	20 seconds
Test strip Storage:	Between +10°C to +32°C
Memory Capacity:	75 Blood Glucose Test Results with Time, Date and Smart Tag of Random, Fasting or Post Prandial (Before or After Meal). Displays the average of Test results for each smart tag.

*Pulsatom*

**Health Care Pvt. Ltd.**

86D2/A, Industrial Suburb, Yeshwantpur, Bangalore 560 022 INDIA  
Tel: +91-80-23579166 Fax: +91-80-23579170 Email: info@pulsatom.com  
Web: www.pulsatom.com